

Digital Sphygmomanometer



DS-300

The symbols on the LCD display

1.Systolic blood pressure

2. Diastolic blood pressure

3.Memory symbol

4. Heartbeat symbol (Flashes during measurement)

5.Average value symbol

7.Date/Time display

9. Irregular heartbeat symbol

11.USER B

13.WHO Function symbol

15.Cuff self-checking function

4. Irregular heartbeat checking

6.Battery low symbol

8. Pulse display

10.Movement error symbol

12. USER A

14.Mute symbol



1.Talking function

2. Double users: 2 x 120 sets memory 3.Cuff self-checking function

5. Average value function

6. Low battery display

7. WHO function

8. Auto power-off

9. External power adapter support 10. Volume adjustment

11. Date/time display

Note: Arm circumference should be measured with a measuring tape in the middle of the relaxed upper arm. Do not force cuff connection into the opening. Make sure the cuff connection is not pushed into the AC adapter port.

Technical Specifications

Wight: 239.3g (Batteries and AC adapter are not included)

Display: 48*65mm [2.56"x1.89"] LCD Digital Display

Size: 118 (W) x 110 (L) x 57 (H) mm 【4.65"(W)x4.33"(L)x2.24"(H)】

Accessories: 1×Main Device, 1×Cuff, 1×Users manual, 1×Warranty card

Operating Conditions: Temperature: 5°C to 40°C;Humidity: 15½ to 93½ RH;

Storage And Shipping Conditions:Temperature: -25°C to 70°C; Humidity:≤ 93′/. RH;

Atmospheric pressure range: 70kPa~106kPa

Measuring method: Oscillometric

Pressure sensor: Resistive Measuring range: 0-280mmHg Pulse: 40 to 199 per minute

Cuff pressure display range:<300mmHg

Memory: Automatically stores the last 120 measurements for 2 users (total 240)

Measuring resolution: 1 mmHg

Accuracy: Pressure within @ 3 mmHg / pulse @ 5 1/2 of the reading

Power source: a) 4*AAA batteries, 1.5 V

b) AC adapter INPUT: 100-240VAC 50/60HZ OUTPUT: 5V DC 1A

Accessories: Wide range rigid cuff 8.7" – 15.7" (22 - 40 cm)

Automatically power off: 60 seconds

Users: Adult

Expected service life of the device and accessories: 5 years

Technical alterations reserved!

Address: 3 Shenton Way #16-08 Shenton House Singapore 068805

Email: international@omega-h-c.com WWW.omega-h-c.com